




# STENTYS

SIMPLE STENT SOLUTIONS

## APPOSITION II TRIAL RESULTS

**R**andomized **C**om**P**arison between the STENTYS Self-**E**x**P**anding **C**oronary **S**tent and a Balloon-Expandable Stent **I**n **A**cu**T**e Myocard**I**al Infarcti**O**N

**DESIGN:** International, prospective, randomized, two-arm, multi-centre trial



### 3 DAY OCT RESULTS

#### % MALAPPOSED STRUTS

STENTYS® Stent	0.51%	p < 0.001
VISION / Driver Stent	5.33%	


#### % PATIENTS WITH STENT MALAPPOSITION\*

STENTYS® Stent	0%	p < 0.001
VISION / Driver Stent	28%	

\* Stent Malapposition defined as more than 5% of struts malapposed under OCT (P. Barlis et al., Eur Heart J (2010) 31 (2) : 165-176)



Perfect Apposition  
of the STENTYS® Stent



Malapposed balloon-expandable stent

#### CONCLUSIONS:


- No STENTYS patients had malapposed stents, whereas 28% of conventional stents were malapposed at 3 days.
- Safety was demonstrated with low MACE at 6 months.
- The STENTYS® Self-Apposing Stent accommodates early changes in AMI vessel anatomy that are caused by thrombus dissolution and vasodilation.

## ACUTE RESULTS<sup>1</sup>


### Primary Endpoint

Stent strut malapposition at 3 days measured by OCT.

Stent Strut Malapposition at 3 days



Patients with Malapposed Stents at 3 days



	STENTYS N = 43	Control N = 37	P Value
Device success	97.6%	100%	NS
Diameter largest balloon used (mm)	2.78	3.32	0.014
At least one balloon inflation	93%	100%	NS
Mean Lumen Area by OCT (mm <sup>2</sup> )			
Post-procedure	7.88 ± 2.32	8.92 ± 2.22	NS
3 days	8.99 ± 2.39	8.81 ± 2.18	NS
3 day in-stent lumen loss by QCA (mm)	-0.11 ± 0.29	0.04 ± 0.21	0.01

**2.3% Clinically-Driven TLR at 6 months<sup>1</sup>**

1. Presented at TCT 2010 by S. Verheye and at PCR 2011 by R.J. van Geuns